

iPSC-Derived Monocytes Generate Functional M1 and M2 Macrophages with Enhanced Cytokine Secretion and Tumor Cell-Killing Activity

lan Hay,¹ Christopher B. Rohde,¹ Matthew Angel¹

¹Factor Bioscience Inc., Cambridge, MA



Introduction and Study Goal



- We have previously demonstrated iPSCs can be generated using an mRNA reprogramming process (Harris *et al.*, Mol Ther, Vol 28 No 4S1, 2020)
- Monocytes and monocyte-derived macrophages possess innate functionality in tumor development; sourcing cells from PBMCs has proven a challenge in research and clinical development
- The goal of this study was to investigate these iPSCs' potential to differentiate into monocytes and monocyte-derived macrophages



© Factor Bioscience Inc.

Macrophages Evident Role in Various Cancers

CD68+ CD163+



Breast Carcinoma

Mwafy, S.E., et al. *J Egypt Natl Canc Inst* **32**, 6 (2020).

FACTOR BIOSCIENCE Colorectal Carcinoma Shabo, I., et al. *Cancer Microenvionment* **7** 61-69 (2014) Pancreatic Carcinoma Gardian, K., et al. J Cancer **3** 285-291 (2012)

CD68+





CD163+

Macrophages

iPSC-Derived Monocytes in 14 days





Day 14 - CD14+ cells in suspension

Day 0 – iPSC colony



- 14-day differentiation protocol yielded CD14+ cells in suspension
 - with large nucleus and round morphology

Yield of CD14+ Monocytes



- 14-day differentiation protocol
- Thawed and plated 10⁶ iPSCs
- Seeded single cell suspension to grow small homogeneous colonies
- Harvest window of ~2 weeks
- Average yield per cm²: 4.1x10⁴ cells
- Average yield per iPSC: 7.1 cells

 iPSCs generated CD14+ monocytes within 18 days from frozen to harvest



	Day 0	Day 14
Tra-1-60	97%	4.1%
Tra-1-81	88%	0.4%

Flow Cytometry of Pluripotency Markers

- Harvested suspension cells negative for pluripotent markers Tra-1-60 and Tra-1-81
- iPSC-Monocytes displayed similar expression to CD14+ PBMCs of key monocyte markers: CD11b, CD13, CD33, CD45

	PBMC- Monocytes	iPSC- Monocytes
CD11b	81%	93%
CD13	99%	99%
CD33	99%	82%
CD45	99%	98%
CD68	0.0%	6.5%
CD80	0.0%	99%
CD206	0.8%	93%
SIRPa	99%	62%

Flow Cytometry Characterization of CD14+ Populations



\bigtriangledown
Ň
U

	PBMC-Monocytes	iPSC-Monocytes
Live/Dead	99%/1.1%	94%/5.8%
CD14	83%	79%
CD45	99%	95%
CD163	13%	3.0%

Flow Cytometry of PBMC- and iPSC-monocytes directly after thawing vials of 5x10⁶ cells/mL in 20% DMSO / 80% RPMI1640 frozen for 1 month

• iPSC-Monocytes capable of cryopreservation with similar results compared to PBMC-Monocytes



iPSC-Monocyte-Derived Macrophages





Day 0

HACIOR bioscience Day 4

- Over 4 days with 50 ng/mL M-CSF, >90% of CD14+ cells
 - adhered, and morphology resembled macrophage shape with prominent vacuoles

X

		Fresh iPSC-Monocytes	Harvest/well	Yield
Polarization	Stimulant	M0 Stimulation	4.17x10 ⁵	83.4%
M1	50 ng/mL IFNy + 10 ng/mL LPS	M1 Stimulation	5.04x10 ⁵	100.8%
M2a	10 ng/mL IL-4	M2a Stimulation	5.43x10 ⁵	108.6%

Cryopreserved Monocytes	Harvest/well	Yield
iPSC-Monocytes (M0)	0.94x10 ⁵	18.8%
PBMC-Monocytes (M0)	0.56x10 ⁵	11.2%

Macrophage yields from plating 500k CD14+ cells per well with 50 ng/mL of M-CSF in a 24-well plate for 4 days

- iPSC-monocytes capable of differentiating into macrophages with high yield
- iPSC-monocytes recover from cryopreservation and differentiate into macrophages with yields comparable to PBMC-monocytes





- iPSC-Macrophages capable of enhanced cytokine secretion after
 - stimulation into both pro- and anti-inflammatory phenotypes











- First reported differentiation and characterization of mRNA-reprogrammed iPSCs into monocyte-derived macrophages
- iPSC-monocytes capable of cryopreservation and recovery with similar yields compared to PBMC alternatives
- iPSC-monocytes able to differentiate into macrophages and polarize in presence of stimulants; cytotoxic towards U2OS cancer cells
- mRNA iPSCs show promise as a perpetual source of therapeutic macrophages





- Optimization of monocyte cryopreservation and recovery
- Characterization of gene expression in polarized macrophages
- Functional assessment of phagocytosis, antigen presentation (T cell activation), tumor infiltration
- Scaling up differentiation process to support pre-clinical and clinical development
- Investigating gene editing for targeting molecules (e.g. chimeric antigen receptors), cytokine expression cassettes, etc.





I.H., C.R. and M.A. are employees of Factor Bioscience Inc.

I.H., C.R. and M.A. are inventors on patents assigned to Factor Bioscience Inc.





Thank you!

ian.hay@factorbio.com